

## Technovent New Product Information

### Vistafix 3 Standard Magnacaps

Standard mini and maxi magnacaps to fit the abutment fixture of the new Vistafix 3 implant system. Use in conjunction with S-Range magnets.

Order Code	Description	Height	Diameter
VF-MC1-S	Mini Magnacap	2.1mm	4.4mm
VF-MC2-S	Maxi Magnacap	2.1mm	5.1mm



### New Auricular Magnacaps for Vistafix Generation 2 and Vistafix 3 Implant Systems

New Auricular maxi magnacaps to fit on the abutment fixture of both Vistafix 2 and Vistafix 3 implant abutments. Must be used in conjunction with auricular magnet (MLL2-OR-S).

Order Code	Description	Height	Diameter
BMC2-OR-S	Maxi Magnacap For Vistafix 2	2.7mm	5.26mm
VF-MC2-OR-S	Maxi Magnacap For Vistafix 3	2.7mm	5.26mm



### New Auricular Magnet

New Auricular magnet designed for the retention of primarily auricular prostheses. S-Range magnetic strength coupled with a mechanical backup reduces lateral displacement. Use in conjunction with either BMC2-OR-S or VF-MC2-OR-S.

Order Code	Description	Height	Diameter
MLL2-OR-S	Auricular Magnet	2.7mm	Fitting surface 5.5mm External at widest point 7.9mm

