

PCSK-10-2026 Platinum Silicone Kit

Silicones

- (1) VST-50: Facial primary silicone
- (1) VST-30C Catalyst only - we simply add a crosslinker quick to mold for soft appliances
- (1) A-RTV-3040 higher durometer, Lining & for fingers, digits, nasal appliance
- (1) A-RTV-3045C Catalyst only - (shore A) 45 durometer

- (1) TS-564 4oz. Sealant: dispersion primary initial seal, after adding extrinsic
- (1) A-564-2 Acetoxy Silicone 2oz: Sealant, paint very thin coat to laminate seal
- (1) MD-564 4oz. matting Agent: remove the gloss after A-564

Impression Materials

- (3) M-517 Soft silicone impression material
- (1) DG-1-S Cartridge gun for dispensing M-517
- (1 SAMPLE) FX-300 impression lube, before impression release and balance ph of skin

Colorants

- (1) FI-200 1 oz white
- (1) FI-202.1 1 oz yellow
- (1) FI-235 1 oz blue
- (1) FI-204 1 oz red
- (1) FI- SK01 1 oz naturelle
- (1) FE-1999 FE color kit
- (1) H-101R Red Flocking
- (1) H-107-P Plum Flocking
- (1) H-108-T Tan Flocking
- (1) H-109-P Pink Flocking

Skin Care

- (1) B-200-R Regular strength Water-based adhesive
- (1) B-200-30 Hydrobond Water based adhesive mild
- (1) B-400 1 oz Silicone adhesive Silicone based adhesive
- (2 SAMPLE WIPES) B-100 Skin Prep Wipes Place on skin prior to adhesive for protection
- (1) B-200-12 Antimicrobial soap Cleansing the skin
- (1) B-200-09 Antimicrobial spray Prevent buildup of Candida on appliance
- (4) B-608-1 oz ProBond adhesive remover Remover for silicone-based adhesive
- (2 SAMPLE WIPES) B-608-M SAMPLE ProBond adhesive remover mitts remover for silicone-based adhesives

PCSK-10 Platinum Silicone Kit (Complete)

Additives

- (1) FE-100 4oz
- (1) A-300-1 Thixotropic additive
- (1) A-330 Gold Primer and 1 Pkg of 5 primer tips

Miscellaneous

Mixing containers

- (12) D-201-P
- (6) PK-6
- (6) PK-8
- (6) PK-16

- (1) 9014 Mixing pad, pkg of 2 (Large)
- (1) G-490-00 Red sable brush
- (1) G-490-0
- (1) G-490-3 “
- (1) G-490-5 “
- (1) G-120 Sculpting Wax ½ sheet (1 sample size)
- (1) G-104 Clay 1 small block (1/4 of the block)

- (1) 5125 polishing wheel kit
- (1) Al-605 Plastic spatula
- (1 PACK, 5 OF EACH SIZE) Sample pack acid brush
- F-901-4

Instructions

- J-914 Textbook “Fundamentals of Facial Prosthetic”
- J-912 Dictionary of Prosthetic Rehabilitation”